


## MICRONIC TUBE SELECTOR

**Manual Tube Selector Tool**  
Rapidly and securely picking of sample storage tubes!



### Ordering information

Catalog number	Description	Unit
MP54007	Micronic Tube Selector Tool	5 pieces

### Superior features

- Easy-to-use tool to manually grab individual tubes from a tube rack without the necessity of lifting them
- Optimized grippers for picking tubes with Screw Caps or Push Caps
- Single-handed operation
- Manual thumb operated ejection
- Convenient grip design
- Minimizing cross-contamination
- Compatible with all Micronic 96-well format tubes
- Wipe clean surface



30 years of experience  
in traceable sample  
storage

Class 7 clean room  
injection molding  
and assembly

Production in  
the Netherlands and  
the United States