



SINGLE CHANNEL PIPETTES

ecopipette Innovation

- THIS PIPETTE IS BUILT TO LAST. MAKING NO COMPROMISE ON QUALITY IN PRODUCT DEVELOPMENT, RAW MATERIALS USED AND MANUFACTURING PROCESSES WE DEVELOPED THE GREENEST PIPETTE AVAILABLE ON THE MARKET. THE ECOPIPETTE IS NOT ONLY PACKAGED USING RECYCLABLE, BIODEGRADABLE MATERIALS, BUT CONSTRUCTED FROM SOME OF THE MOST RENEWABLE DURABLE RESOURCES AVAILABLE TO ENSURE AN ADVANCED LIFE SPAN AND GREENEST AFTER SERVICE RECYCLING.
- CAPP ECOTRADE™ PROGRAM: RETURN ANY BRAND OF OLD PIPETTES IN EXCHANGE FOR DISCOUNTED NEW ECOPIPETTES™. AT CAPP APS, WE WILL ENSURE AN ENVIRONMENTALLY SAFE WASTE DISCHARGE OF ALL OLD PIPETTES. ASK YOUR FAVOURITE AUTHORIZED CAPP DISTRIBUTOR FOR MORE INFORMATION ABOUT OUR ECOTRADE™ PROGRAMME.

CAPP ECOPIPETTE™ UNIQUE ADVANTAGES:

- LOCKABLE VOLUME CONTROLLER REQUIRING 50% FEWER TURNS
- DISPLAY IS FULLY VISIBLE DURING HAND OPERATIONS
- REDUCED FORCE REQUIREMENT FOR ASPIRATION AND BLOW-OUT
- THE ONLY PIPETTE WITH FIXED AND VARIABLE EXCHANGEABLE COUNTER MECHANISM
- LARGE PAD FOR EASE OF TIP EJECTION
- FULLY AUTOCLAVABLE PIPETTE BODY (NO DISASSEMBLY REQUIRED)
- DOUBLE O-RING FOR OPTIMAL TIP SEAL
- COMPATIBLE WITH ALL MAJOR PIPETTE TIP BRANDS
- ROBUST DESIGN INCLUDES METAL TIP BASE
- LOWEST INCIDENCE OF BREAKAGE RESULTING FROM ACCIDENTAL DROPPING OR TIPS MOUNTING.





Ordering Info

ecopipette variable volume

CAT. NO.	DESCRIPTION	INACCURACY %	IMPRECISION %	COLOR
C02-1S	CAPP PIPETTE, VARIABLE VOL. 0.2-2 UL, ECOPIPETTE	12.00 / 2.00	6.00 / 1.00	
C10-1	CAPP PIPETTE, VARIABLE VOL. 0.5-10 UL, ECOPIPETTE	2.50 / 1.00	1.80 / 0.50	
C20-1	CAPP PIPETTE, VARIABLE VOL. 2-20 UL, ECOPIPETTE	5.00 / 1.00	1.50 / 0.30	
C50-1	CAPP PIPETTE, VARIABLE VOL. 5-50 UL, ECOPIPETTE	2.00 / 1.00	1.20 / 0.30	
C100-1	CAPP PIPETTE, VARIABLE VOL. 10-100 UL, ECOPIPETTE	2.50 / 0.70	1.00 / 0.20	
C200-1	CAPP PIPETTE, VARIABLE VOL. 20-200 UL, ECOPIPETTE	2.50 / 0.60	0.70 / 0.20	
C1000-1	CAPP PIPETTE, VARIABLE VOL. 100-1000 UL, ECOPIPETTE	2.00 / 0.80	0.60 / 0.20	
C5000-1	CAPP PIPETTE, VARIABLE VOL. 1-5 ML, ECOPIPETTE	1.20 / 0.60	0.30 / 0.15	

ecopipette fixed volume (customized volumes available on request)

CAT. NO.	DESCRIPTION	INACCURACY %	IMPRECISION %
C02-1S-001	CAPP PIPETTE, FIXED VOLUME 1 UL , ECOPIPETTE	2.50	1.80
C10-1-005	CAPP PIPETTE, FIXED VOLUME 5 UL, ECOPIPETTE	1.20	0.80
C10-1-010	CAPP PIPETTE, FIXED VOLUME 10 UL, ECOPIPETTE	1.00	0.50
C50-1-020	CAPP PIPETTE, FIXED VOLUME 20 UL, ECOPIPETTE	1.00	0.30
C50-1-025	CAPP PIPETTE, FIXED VOLUME 25 UL, ECOPIPETTE	0.80	0.35
C100-1-030	CAPP PIPETTE, FIXED VOLUME 30 UL, ECOPIPETTE	1.16	0.40
C100-1-050	CAPP PIPETTE, FIXED VOLUME 50 UL, ECOPIPETTE	1.00	0.30
C200-1-100	CAPP PIPETTE, FIXED VOLUME 100 UL , ECOPIPETTE	0.70	0.20
C200-1-200	CAPP PIPETTE, FIXED VOLUME 200 UL, ECOPIPETTE	0.60	0.20
C1000-1-500	CAPP PIPETTE, FIXED VOLUME 500 UL, ECOPIPETTE	0.80	0.15
C1000-1-1000	CAPP PIPETTE, FIXED VOLUME 1000 UL, ECOPIPETTE	0.80	0.20
C5000-1A-5000	CAPP PIPETTE, FIXED VOLUME 5000 UL, COMFORT	0.60	0.15

CXXX-1-YYY "OTHER FIXED VOL. PIPETTES(NOTE: XXX=MAX PIPETTE VOLUME; YYY=DESIRED VOLUME)"

ecopipette Starter Packs

CAT. No.	DESCRIPTION
SKE-03-XX	THREE PIPETTES FREE OF CHOICE, STAND, RACKS OF TIPS AND TIMER